

Antibodies Product Datasheet

Product Name:	CD22 Antibody
Catalog Number:	B8230
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BL-CAM. AA range:776-825
Concentration:	1 mg/ml
Gene Name:	CD22
Protein Name:	B-cell receptor CD22
Other Name:	CD22; SIGLEC2; B-cell receptor CD22; B-lymphocyte cell adhesion molecule; BL-CAM; Sialic acid-binding Ig-like lectin 2; Siglec-2; T-cell surface antigen Leu-14; CD antigen CD22
Human Swiss Prot No.:	P20273
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD22 Polyclonal Antibody detects endogenous levels of CD22 protein.
Storage and Stability:	-20°C/1 year