

Antibodies Product Datasheet

Product Name:	Fer Antibody
Catalog Number:	B8100
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FER. AA range:371-420
Concentration:	1 mg/ml
Gene Name:	FER
Protein Name:	Tyrosine-protein kinase Fer
Other Name:	FER; TYK3; Tyrosine-protein kinase Fer; Feline encephalitis virus-related kinase FER; Fujinami poultry sarcoma/Feline sarcoma-related protein Fer; Proto-oncogene c-Fer; Tyrosine kinase 3; p94-Fer
Human Swiss Prot No.:	P16591
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fer Polyclonal Antibody detects endogenous levels of Fer protein.
Storage and Stability:	-20°C/1 year