

## **Assay Biotechnology Company**

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

## **Antibodies Product Datasheet**

**Product Name:** Stat4 Antibody

Catalog Number: B7225

Reactivity: Human, Mouse, Rat

**Applications:** WB, IHC-p, IF(paraffin section), IP, ELISA

**Dilutions:** Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/5000. Not yet tested in

other applications.

**Modification:** Unmodified/Total

**Source:** Polyclonal Rabbit

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Immunogen:** The antiserum was produced against synthesized peptide derived from

human STAT4. AA range:660-709

Concentration: 1 mg/ml

Gene Name: STAT4

**Protein Name:** Signal transducer and activator of transcription 4

Other Name: STAT4; Signal transducer and activator of transcription 4

Human Swiss Prot No.: Q14765

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

**Specificity:** Stat4 Polyclonal Antibody detects endogenous levels of Stat4 protein.

Storage and Stability: -20°C/1 year