

Antibodies Product Datasheet

Product Name:	Jun B Antibody
Catalog Number:	B7135
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), IP, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human JunB. AA range:226-275
Concentration:	1 mg/ml
Gene Name:	JUNB
Protein Name:	Transcription factor jun-B
Other Name:	JUNB; Transcription factor jun-B
Human Swiss Prot No.:	P17275
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Jun B Polyclonal Antibody detects endogenous levels of Jun B protein.
Storage and Stability:	-20°C/1 year