

Antibodies Product Datasheet

Product Name:	JAK2 Antibody
Catalog Number:	B7126
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human JAK2. AA range:981-1030
Concentration:	1 mg/ml
Gene Name:	JAK2
Protein Name:	Tyrosine-protein kinase JAK2
Other Name:	JAK2; Tyrosine-protein kinase JAK2; Janus kinase 2; JAK-2
Human Swiss Prot No.:	O60674
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	JAK2 Polyclonal Antibody detects endogenous levels of JAK2 protein.
Storage and Stability:	-20°C/1 year