

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: GluR-1 Antibody

Catalog Number: B7094

Reactivity: Human, Mouse, Rat

Applications: WB, IHC-p, IF/ICC, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in

other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human GluR1. AA range:816-865

Concentration: 1 mg/ml

Gene Name: GRIA1

Protein Name: Glutamate receptor 1

Other Name: GRIA1; GLUH1; GLUR1; Glutamate receptor 1; GluR-1; AMPA-selective

glutamate receptor 1; GluR-A; GluR-K1; Glutamate receptor ionotropic;

AMPA 1; GluA1

Human Swiss Prot No.: P42261

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: GluR-1 Polyclonal Antibody detects endogenous levels of GluR-1

protein.

Storage and Stability: -20°C/1 year