

## Antibodies Product Datasheet

<b>Product Name:</b>	ATF-2 Antibody
<b>Catalog Number:</b>	B7014
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), IP, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human ATF2. AA range:35-84
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ATF2
<b>Protein Name:</b>	Cyclic AMP-dependent transcription factor ATF-2
<b>Other Name:</b>	ATF2; CREB2; CREBP1; Cyclic AMP-dependent transcription factor ATF-2; cAMP-dependent transcription factor ATF-2; Activating transcription factor 2; Cyclic AMP-responsive element-binding protein 2; CREB-2; cAMP-responsive element-binding protein
<b>Human Swiss Prot No.:</b>	P15336
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	ATF-2 Polyclonal Antibody detects endogenous levels of ATF-2 protein.
<b>Storage and Stability:</b>	-20°C/1 year