

Antibodies Product Datasheet

Product Name:	AR Antibody
Catalog Number:	B7009
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Androgen Receptor. AA range:621-670
Concentration:	1 mg/ml
Gene Name:	AR
Protein Name:	Androgen receptor
Other Name:	AR; DHTR; NR3C4; Androgen receptor; Dihydrotestosterone receptor; Nuclear receptor subfamily 3 group C member 4
Human Swiss Prot No.:	P10275
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AR Polyclonal Antibody detects endogenous levels of AR protein.
Storage and Stability:	-20°C/1 year