

Antibodies Product Datasheet

Product Name:	Pdcd-4 Antibody
Catalog Number:	B1174
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PDCD4. AA range:33-82
Concentration:	1 mg/ml
Gene Name:	PDCD4
Protein Name:	Programmed cell death protein 4
Other Name:	PDCD4; H731; Programmed cell death protein 4; Neoplastic transformation inhibitor protein; Nuclear antigen H731-like; Protein 197/15a
Human Swiss Prot No.:	Q53EL6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Pdcd-4 Polyclonal Antibody detects endogenous levels of Pdcd-4 protein.
Storage and Stability:	-20°C/1 year