

Antibodies Product Datasheet

Product Name:	AR- β 2 Antibody
Catalog Number:	B1131
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Adrenergic Receptor beta2. AA range:321-370
Concentration:	1 mg/ml
Gene Name:	ADRB2
Protein Name:	Beta-2 adrenergic receptor
Other Name:	ADRB2; ADRB2R; B2AR; Beta-2 adrenergic receptor; Beta-2 adrenoceptor; Beta-2 adrenoceptor
Human Swiss Prot No.:	P07550
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AR- β 2 Polyclonal Antibody detects endogenous levels of AR- β 2 protein.
Storage and Stability:	-20°C/1 year