

Antibodies Product Datasheet

Product Name:	FAS Antibody
Catalog Number:	B0927
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	WB 1:500-2000, IF 1:50-300, IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FAS. AA range:257-306
Concentration:	1 mg/ml
Gene Name:	FAS
Protein Name:	Tumor necrosis factor receptor superfamily member 6
Other Name:	FAS; APT1; FAS1; TNFRSF6; Tumor necrosis factor receptor superfamily member 6; Apo-1 antigen; Apoptosis-mediating surface antigen FAS; FASLG receptor; CD antigen CD95
Human Swiss Prot No.:	P25445
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FAS Polyclonal Antibody detects endogenous levels of FAS protein.
Storage and Stability:	-20°C/1 year