

Antibodies Product Datasheet

Product Name:	Met Antibody
Catalog Number:	B0871
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MET. AA range:976-1025
Concentration:	1 mg/ml
Gene Name:	MET
Protein Name:	Hepatocyte growth factor receptor
Other Name:	MET; Hepatocyte growth factor receptor; HGF receptor; HGF/SF receptor; Proto-oncogene c-Met; Scatter factor receptor; SF receptor; Tyrosine-protein kinase Met
Human Swiss Prot No.:	P08581
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Met Polyclonal Antibody detects endogenous levels of Met protein.
Storage and Stability:	-20°C/1 year