

Antibodies Product Datasheet

Product Name:	CD71/TfR Antibody
Catalog Number:	B0852
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD71/TfR. AA range:15-64
Concentration:	1 mg/ml
Gene Name:	TFRC
Protein Name:	Transferrin receptor protein 1
Other Name:	TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD antigen CD71
Human Swiss Prot No.:	P02786
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD71 Polyclonal Antibody detects endogenous levels of CD71 protein.
Storage and Stability:	-20°C/1 year