

Antibodies Product Datasheet

Product Name:	Caspase-3 Antibody
Catalog Number:	B0835
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Caspase 3. AA range:116-165
Concentration:	1 mg/ml
Gene Name:	CASP3
Protein Name:	Caspase3
Other Name:	CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1
Human Swiss Prot No.:	P42574
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Caspase-3 Polyclonal Antibody detects endogenous levels of Caspase-3 protein.
Storage and Stability:	-20°C/1 year