

Antibodies Product Datasheet

Product Name:	Akt Antibody
Catalog Number:	B0817
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Akt. AA range:292-341
Concentration:	1 mg/ml
Gene Name:	AKT1/AKT2/AKT3
Protein Name:	RAC-alpha serine/threonine-protein kinase
Other Name:	AKT1; PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha
Human Swiss Prot No.:	P31749
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Akt Polyclonal Antibody detects endogenous levels of Akt protein.
Storage and Stability:	-20°C/1 year