
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Raf-B Antibody |
| Catalog Number: | B0781 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-p, IF(paraffin section), ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human B-RAF. AA range:576-625 |
| Concentration: | 1 mg/ml |
| Gene Name: | BRAF |
| Protein Name: | Serine/threonine-protein kinase B-raf |
| Other Name: | BRAF; BRAF1; RAFB1; Serine/threonine-protein kinase B-raf; Proto-oncogene B-Raf; p94; v-Raf murine sarcoma viral oncogene homolog B1 |
| Human Swiss Prot No.: | P15056 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Raf-B Polyclonal Antibody detects endogenous levels of Raf-B protein. |
| Storage and Stability: | -20°C/1 year |