

## Antibodies Product Datasheet

<b>Product Name:</b>	SHIP-1 Antibody
<b>Catalog Number:</b>	B0738
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human SHIP1. AA range:987-1036
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	INPP5D
<b>Protein Name:</b>	Phosphatidylinositol 3,4,5-trisphosphate 5-phosphatase 1
<b>Other Name:</b>	INPP5D; SHIP; SHIP1; Phosphatidylinositol 3; 4,5-trisphosphate 5-phosphatase 1; Inositol polyphosphate-5-phosphatase of 145 kDa; SIP-145; SH2 domain-containing inositol 5'-phosphatase 1; SH2 domain-containing inositol phosphatase 1; SHIP-1;
<b>Human Swiss Prot No.:</b>	Q92835
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	SHIP-1 Polyclonal Antibody detects endogenous levels of SHIP-1 protein.
<b>Storage and Stability:</b>	-20°C/1 year