

Antibodies Product Datasheet

Product Name:	Fyn Antibody
Catalog Number:	B0430
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Fyn. AA range:488-537
Concentration:	1 mg/ml
Gene Name:	FYN
Protein Name:	Tyrosine-protein kinase Fyn
Other Name:	FYN; Tyrosine-protein kinase Fyn; Proto-oncogene Syn; Proto-oncogene c-Fyn; Src-like kinase; SLK; p59-Fyn
Human Swiss Prot No.:	P06241
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fyn Polyclonal Antibody detects endogenous levels of Fyn protein.
Storage and Stability:	-20°C/1 year