
Antibodies Product Datasheet

Product Name:	Daxx Antibody
Catalog Number:	B0419
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Daxx. AA range:634-683
Concentration:	1 mg/ml
Gene Name:	DAXX
Protein Name:	Death domain-associated protein 6
Other Name:	DAXX; BING2; DAP6; Death domain-associated protein 6; Daxx; hDaxx; ETS1-associated protein 1; EAP1; Fas death domain-associated protein Q9UER7
Human Swiss Prot No.:	Q9UER7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Daxx Polyclonal Antibody detects endogenous levels of Daxx protein.
Storage and Stability:	-20°C/1 year