

## Antibodies Product Datasheet

<b>Product Name:</b>	4E-BP1 Antibody
<b>Catalog Number:</b>	B0401
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human 4E-BP1. AA range:30-79
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	EIF4EBP1
<b>Protein Name:</b>	Eukaryotic translation initiation factor 4E-binding protein 1
<b>Other Name:</b>	EIF4EBP1; Eukaryotic translation initiation factor 4E-binding protein 1; 4E-BP1; eIF4E-binding protein 1; Phosphorylated heat- and acid-stable protein regulated by insulin 1; PHAS-I
<b>Human Swiss Prot No.:</b>	Q13541
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	4E-BP1 Polyclonal Antibody detects endogenous levels of 4E-BP1 protein.
<b>Storage and Stability:</b>	-20°C/1 year