

## Antibodies Product Datasheet

<b>Product Name:</b>	Caspase-8 Antibody
<b>Catalog Number:</b>	B0059
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IF/ICC, WB, IHC-p, ELISA
<b>Dilutions:</b>	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Caspase 8. AA range:313-362
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CASP8
<b>Protein Name:</b>	Caspase8
<b>Other Name:</b>	CASP8; MCH5; Caspase-8; CASP-8; Apoptotic cysteine protease; Apoptotic protease Mch-5; CAP4; FADD-homologous ICE/ced-3-like protease; FADD-like ICE; FLICE; ICE-like apoptotic protease 5; MORT1-associated ced-3 homolog; MACH
<b>Human Swiss Prot No.:</b>	Q14790
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Caspase-8 Polyclonal Antibody detects endogenous levels of Caspase-8 protein.
<b>Storage and Stability:</b>	-20°C/1 year