

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Smad2 Antibody |
| Catalog Number: | B0030 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-p, IF/ICC, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human Smad2. AA range:418-467 |
| Concentration: | 1 mg/ml |
| Gene Name: | SMAD2 |
| Protein Name: | Mothers against decapentaplegic homolog 2 |
| Other Name: | SMAD2; MADH2; MADR2; Mothers against decapentaplegic homolog 2; MAD homolog 2; Mothers against DPP homolog 2; JV18-1; Mad-related protein 2; hMAD-2; SMAD family member 2; SMAD 2; Smad2; hSMAD2 |
| Human Swiss Prot No.: | Q15796 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Smad2 Polyclonal Antibody detects endogenous levels of Smad2 protein. |
| Storage and Stability: | -20°C/1 year |