

Antibodies Product Datasheet

Product Name:	HSP90 β Antibody
Catalog Number:	B0013
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IP, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HSP90B. AA range:226-275
Concentration:	1 mg/ml
Gene Name:	HSP90AB1
Protein Name:	Heat shock protein HSP 90-beta
Other Name:	HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
Human Swiss Prot No.:	P08238
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HSP90 β Polyclonal Antibody detects endogenous levels of HSP90 β protein.
Storage and Stability:	-20°C/1 year