

Antibodies Product Datasheet

Product Name:	MCM3 (Phospho-Ser112) Antibody
Catalog Number:	A9418
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human MCM3 (Ser112)
Concentration:	1 mg/ml
Gene Name:	MCM3
Protein Name:	MCM3 (Ser112)
Other Name:	DNA replication licensing factor MCM3 (EC 3.6.4.12) (DNA polymerase alpha holoenzyme-associated protein P1) (P1-MCM3) (RLF subunit beta) (p102)
Human Swiss Prot No.:	P25205
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human MCM3 (phospho-Ser112)
Storage and Stability:	-20°C/1 year