
Antibodies Product Datasheet

Product Name:	INCENP (Phospho-Thr59) Antibody
Catalog Number:	A9093
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human INCENP (Phospho-Thr59)
Concentration:	1 mg/ml
Gene Name:	INCENP
Protein Name:	INCENP (Phospho-Thr59)
Other Name:	Inner centromere protein
Human Swiss Prot No.:	Q9NQS7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.186% sodium azide.
Specificity:	This antibody detects endogenous levels of INCENP (Phospho-Thr59) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year