

Antibodies Product Datasheet

Product Name:	TGFBR1 (Phospho-Thr204) Antibody
Catalog Number:	A9060
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human TGFBR1 (Phospho-Thr204)
Concentration:	1 mg/ml
Gene Name:	TGFBR1 ALK5 SKR4
Protein Name:	TGFBR1 (Phospho-Thr204)
Other Name:	TGF-beta receptor type-1 (TGFR-1) (EC 2.7.11.30) (Activin A receptor type II-like protein kinase of 53kD) (Activin receptor-like kinase 5) (ALK-5) (ALK5) (Serine/threonine-protein kinase receptor R4) (SKR4) (TGF-beta type I receptor) (Transforming growth
Human Swiss Prot No.:	P36897
Formulation:	Liquid in PBS containing 50% glycerol, and 0.103% sodium azide.
Specificity:	This antibody detects endogenous levels of TGFBR1 (Phospho-Thr204) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year