
Antibodies Product Datasheet

Product Name:	GRK1/2 (Phospho-Tyr13) Antibody
Catalog Number:	A8508
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GRK1/2 (Phospho-Tyr13)
Concentration:	1 mg/ml
Gene Name:	ADRBK1 BARK BARK1 GRK2
Protein Name:	GRK1/2 (Phospho-Tyr13)
Other Name:	Beta-adrenergic receptor kinase 1 (Beta-ARK-1) (EC 2.7.11.15) (G-protein coupled receptor kinase 2)
Human Swiss Prot No.:	P25098
Formulation:	Liquid in PBS containing 50% glycerol, and 0.216% sodium azide.
Specificity:	This antibody detects endogenous levels of GRK1/2 (Phospho-Tyr13) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year