

Antibodies Product Datasheet

Product Name:	Cyclin D2 (Phospho-Thr280) Antibody
Catalog Number:	A8336
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cyclin D2 around the phosphorylation site of Thr280. AA range:240-289
Concentration:	1 mg/ml
Gene Name:	CCND2
Protein Name:	G1/S-specific cyclin-D2
Other Name:	CCND2; G1/S-specific cyclin-D2
Human Swiss Prot No.:	P30279
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Cyclin D2 (T280) Polyclonal Antibody detects endogenous levels of Cyclin D2 protein only when phosphorylated at T280.
Storage and Stability:	-20°C/1 year