

Antibodies Product Datasheet

Product Name:	IL-4R α Antibody
Catalog Number:	R12-2211
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL-4R α . AA range:205-254
Concentration:	1 mg/ml
Gene Name:	IL4R
Protein Name:	Interleukin-4 receptor subunit alpha
Other Name:	IL4R; IL4RA; 582J2.1; Interleukin-4 receptor subunit alpha; IL-4 receptor subunit alpha; IL-4R subunit alpha; IL-4R-alpha; IL-4RA; CD antigen CD124
Human Swiss Prot No.:	P24394
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-4R α Polyclonal Antibody detects endogenous levels of IL-4R α protein.
Storage and Stability:	-20°C/1 year