

Antibodies Product Datasheet

Product Name:	AIRE-1 (Phospho-Ser156) Antibody
Catalog Number:	A8286
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AIRE around the phosphorylation site of Ser156. AA range:126-175
Concentration:	1 mg/ml
Gene Name:	AIRE
Protein Name:	Autoimmune regulator
Other Name:	AIRE; APECED; Autoimmune regulator; Autoimmune polyendocrinopathy candidiasis ectodermal dystrophy protein; APECED protein
Human Swiss Prot No.:	O43918
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-AIRE-1 (S156) Polyclonal Antibody detects endogenous levels of AIRE-1 protein only when phosphorylated at S156.
Storage and Stability:	-20°C/1 year