

Antibodies Product Datasheet

Product Name:	FAK (Phospho-Ser732) Antibody
Catalog Number:	A8268
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FAK (Phospho-Ser732)
Concentration:	1 mg/ml
Gene Name:	PTK2 FAK FAK1
Protein Name:	FAK (Phospho-Ser732)
Other Name:	Focal adhesion kinase 1 (FADK 1) (EC 2.7.10.2) (Focal adhesion kinase-related nonkinase) (FRNK) (Protein phosphatase 1 regulatory subunit 71) (PPP1R71) (Protein-tyrosine kinase 2) (p125FAK) (pp125FAK)
Human Swiss Prot No.:	Q05397
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous phospho levels of FAK (Phospho-Ser732) at Human:S732, Mouse:S770, Rat:S732
Storage and Stability:	-20°C/1 year