

Antibodies Product Datasheet

Product Name:	IKK γ (Phospho-Ser376) Antibody
Catalog Number:	A8114
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IKK-gamma around the phosphorylation site of Ser376. AA range:342-391
Concentration:	1 mg/ml
Gene Name:	IKBKG
Protein Name:	NF-kappa-B essential modulator
Other Name:	IKBKG; FIP3; NEMO; NF-kappa-B essential modulator; NEMO; FIP-3; I κ B kinase-associated protein 1; IKKAP1; Inhibitor of nuclear factor kappa-B kinase subunit gamma; I-kappa-B kinase subunit gamma; IKK-gamma; IKKG; I κ B kinase subunit gamma; NF
Human Swiss Prot No.:	Q9Y6K9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-IKK γ (S376) Polyclonal Antibody detects endogenous levels of IKK γ protein only when phosphorylated at S376.
Storage and Stability:	-20°C/1 year