

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | EGFR (Phospho-Tyr998) Antibody |
| Catalog Number: | A8026 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | IHC-p, IF(paraffin section), WB |
| Dilutions: | IHC-p 1:50-200, WB 1:500-2000 |
| Modification: | Phospho Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Immunogen: | Synthesized peptide derived from human EGFR (Phospho-Tyr998) |
| Concentration: | 1 mg/ml |
| Gene Name: | EGFR ERBB ERBB1 HER1 |
| Protein Name: | EGFR (Phospho-Tyr998) |
| Other Name: | Epidermal growth factor receptor (EC 2.7.10.1) (Proto-oncogene c-ErbB-1) (Receptor tyrosine-protein kinase erbB-1) |
| Human Swiss Prot No.: | P00533 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous phospho levels of EGFR (Phospho-Tyr998) at Human:Y998, Mouse:Y1000, Rat:Y999 |
| Storage and Stability: | -20°C/1 year |