

## Antibodies Product Datasheet

<b>Product Name:</b>	EGFR (Phospho-Ser1071) Antibody
<b>Catalog Number:</b>	A8025
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human EGFR around the phosphorylation site of Ser1071. AA range:1036-1085
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	EGFR
<b>Protein Name:</b>	Epidermal growth factor receptor
<b>Other Name:</b>	EGFR; ERBB; ERBB1; HER1; Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1
<b>Human Swiss Prot No.:</b>	P00533
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-EGFR (S1071) Polyclonal Antibody detects endogenous levels of EGFR protein only when phosphorylated at S1071.
<b>Storage and Stability:</b>	-20°C/1 year