Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: CREB-1 (Phospho-Thr100) Antibody

Catalog Number: A8021

Reactivity: Human, Mouse, Rat

Applications: WB, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other

applications.

Modification: Phospho Specific

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human CREB around the phosphorylation site of Thr100. AA range:71-

120

Concentration: 1 mg/ml

Gene Name: CREB1

Protein Name: Cyclic AMP-responsive element-binding protein 1

Other Name: CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1;

cAMP-responsive element-binding protein 1

Human Swiss Prot No.: P16220

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Phospho-CREB-1 (T100) Polyclonal Antibody detects endogenous levels

of CREB-1 protein only when phosphorylated at T100.

Storage and Stability: -20°C/1 year

Last updated December 2022