

## Antibodies Product Datasheet

<b>Product Name:</b>	p53 (Phospho-Ser6) Antibody
<b>Catalog Number:</b>	A7185
<b>Reactivity:</b>	Human, Monkey
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), IP, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human p53 around the phosphorylation site of Ser6. AA range:1-50
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	TP53
<b>Protein Name:</b>	Cellular tumor antigen p53
<b>Other Name:</b>	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
<b>Human Swiss Prot No.:</b>	P04637
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-p53 (S6) Polyclonal Antibody detects endogenous levels of p53 protein only when phosphorylated at S6.
<b>Storage and Stability:</b>	-20°C/1 year