

## Antibodies Product Datasheet

<b>Product Name:</b>	Jun B (Phospho-Ser79) Antibody
<b>Catalog Number:</b>	A7136
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC-p, IF(paraffin section), IP, ELISA
<b>Dilutions:</b>	WB 1:500-2000 ,Immunohistochemistry: 1/100 - 1/300. Immunoprecipitation: 2-5 ug/mg lysate. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human JunB around the phosphorylation site of Ser79. AA range:46-95
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	JUNB
<b>Protein Name:</b>	Transcription factor jun-B
<b>Other Name:</b>	JUNB; Transcription factor jun-B
<b>Human Swiss Prot No.:</b>	P17275
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-Jun B (S79) Polyclonal Antibody detects endogenous levels of Jun B protein only when phosphorylated at S79.
<b>Storage and Stability:</b>	-20°C/1 year