

---

## Antibodies Product Datasheet

<b>Product Name:</b>	AP-1 (Phospho-Ser63) Antibody
<b>Catalog Number:</b>	A7045
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IF/ICC, WB, IHC-p, ELISA
<b>Dilutions:</b>	IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human c-Jun around the phosphorylation site of Ser63. AA range:31-80
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	JUN
<b>Protein Name:</b>	Transcription factor AP-1
<b>Other Name:</b>	JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39
<b>Human Swiss Prot No.:</b>	P05412
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-AP-1 (S63) Polyclonal Antibody detects endogenous levels of AP-1 protein only when phosphorylated at S63.
<b>Storage and Stability:</b>	-20°C/1 year