

Antibodies Product Datasheet

Product Name:	Cdc25A (Phospho-Ser124) Antibody
Catalog Number:	A0854
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDC25A around the phosphorylation site of Ser124. AA range:90-139
Concentration:	1 mg/ml
Gene Name:	CDC25A
Protein Name:	M-phase inducer phosphatase 1
Other Name:	CDC25A; M-phase inducer phosphatase 1; Dual specificity phosphatase Cdc25A
Human Swiss Prot No.:	P30304
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Cdc25A (S124) Polyclonal Antibody detects endogenous levels of Cdc25A protein only when phosphorylated at S124.
Storage and Stability:	-20°C/1 year