Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: Akt (Phospho-Tyr326) Antibody

Catalog Number: A0817

Reactivity: Human, Mouse, Rat

Applications: WB, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other

applications.

Modification: Phospho Specific

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human Akt around the phosphorylation site of Tyr326. AA range:292-341

Concentration: 1 mg/ml

Gene Name: AKT1/AKT2/AKT3

Protein Name: RAC-alpha serine/threonine-protein kinase

Other Name: AKT1; PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein

kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-

Akt; RAC-PK-alpha

Human Swiss Prot No.: P31749

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Phospho-Akt (Y326) Polyclonal Antibody detects endogenous levels of

Akt protein only when phosphorylated at Y326.

Storage and Stability: -20°C/1 year