

Antibodies Product Datasheet

Product Name:	Vav (Phospho-Tyr160) Antibody
Catalog Number:	A0593
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human VAV1 around the phosphorylation site of Tyr160. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	VAV1
Protein Name:	Proto-oncogene vav
Other Name:	VAV1; VAV; Proto-oncogene vav
Human Swiss Prot No.:	P15498
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Vav (Y160) Polyclonal Antibody detects endogenous levels of Vav protein only when phosphorylated at Y160.
Storage and Stability:	-20°C/1 year