

## Antibodies Product Datasheet

<b>Product Name:</b>	HSP90 $\beta$ (Phospho-Ser254) Antibody
<b>Catalog Number:</b>	A0013
<b>Reactivity:</b>	Human, Mouse, Rat, Monkey
<b>Applications:</b>	WB, IHC-p, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human HSP90B around the phosphorylation site of Ser254. AA range:226-275
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	HSP90AB1
<b>Protein Name:</b>	Heat shock protein HSP 90-beta
<b>Other Name:</b>	HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
<b>Human Swiss Prot No.:</b>	P08238
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-HSP90 $\beta$ (S254) Polyclonal Antibody detects endogenous levels of HSP90 $\beta$ protein only when phosphorylated at S254.
<b>Storage and Stability:</b>	-20°C/1 year