

Antibodies Product Datasheet

Product Name:	BID (p11, Cleaved-Ser100) Antibody
Catalog Number:	L0476
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BID (p11, Cleaved-Ser100)
Concentration:	1 mg/ml
Gene Name:	BID
Protein Name:	BID (p11, Cleaved-Ser100)
Other Name:	BH3-interacting domain death agonist (p22 BID) (BID) [Cleaved into: BH3-interacting domain death agonist p15 (p15 BID); BH3-interacting domain death agonist p13 (p13 BID); BH3-interacting domain death agonist p11 (p11 BID)]
Human Swiss Prot No.:	P55957
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse BID (p11, Cleaved-Ser100, protein was cleaved amino acid sequence between 99-100)
Storage and Stability:	-20°C/1 year