

Antibodies Product Datasheet

Product Name:	THRB (heavy chain, Cleaved-Ile364) Antibody
Catalog Number:	L0456
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human THRB (heavy chain, Cleaved-Ile364)
Concentration:	1 mg/ml
Gene Name:	F2
Protein Name:	THRB (heavy chain, Cleaved-Ile364)
Other Name:	Prothrombin (EC 3.4.21.5) (Coagulation factor II) [Cleaved into: Activation peptide fragment 1; Activation peptide fragment 2; Thrombin light chain; Thrombin heavy chain]
Human Swiss Prot No.:	P00734
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat THRB (heavy chain, Cleaved-Ile364, protein was cleaved amino acid sequence between 363-364)
Storage and Stability:	-20°C/1 year