

Antibodies Product Datasheet

Product Name:	Cleaved-MPO 89k (A49) Antibody
Catalog Number:	L0124
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Cleaved-MPO 89k (A49) . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	MPO
Protein Name:	Myeloperoxidase
Other Name:	MPO; Myeloperoxidase; MPO
Human Swiss Prot No.:	P05164
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-MPO 89k (A49) Polyclonal Antibody detects endogenous levels of fragment of activated MPO 89k protein resulting from cleavage adjacent to A49.
Storage and Stability:	-20°C/1 year