

Antibodies Product Datasheet

Product Name:	Cleaved-Caspase-6 p18 (D179) Antibody
Catalog Number:	L0106
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Caspase 6. AA range:130-179
Concentration:	1 mg/ml
Gene Name:	CASP6
Protein Name:	Caspase6
Other Name:	CASP6; MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2
Human Swiss Prot No.:	P55212
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Caspase-6 p18 (D179) Polyclonal Antibody detects endogenous levels of fragment of activated Caspase-6 p18 protein resulting from cleavage adjacent to D179.
Storage and Stability:	-20°C/1 year