

Antibodies Product Datasheet

Product Name:	Cleaved-Aggrecan (D369) Antibody
Catalog Number:	L0101
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Aggrecan. AA range:320-369
Concentration:	1 mg/ml
Gene Name:	ACAN
Protein Name:	Aggrecan core protein
Other Name:	ACAN; AGC1; CSPG1; MSK16; Aggrecan core protein; Cartilage- specific proteoglycan core protein; CSPCP; Chondroitin sulfate proteoglycan core protein 1; Chondroitin sulfate proteoglycan 1
Human Swiss Prot No.:	P16112
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Aggrecan (D369) Polyclonal Antibody detects endogenous levels of fragment of activated Aggrecan protein resulting from cleavage adjacent to D369.
Storage and Stability:	-20°C/1 year