
Antibodies Product Datasheet

Product Name:	RUNX1 (Acetyl-Lys24) Antibody
Catalog Number:	D12316
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human RUNX1 (Acetyl Lys24)
Concentration:	1 mg/ml
Gene Name:	RUNX1 AML1 CBFA2
Protein Name:	RUNX1 (Acetyl Lys24)
Other Name:	Runt-related transcription factor 1 (Acute myeloid leukemia 1 protein) (Core-binding factor subunit alpha-2) (CBF-alpha-2) (Oncogene AML-1) (Polyomavirus enhancer-binding protein 2 alpha B subunit) (PEA2-alpha B) (PEBP2-alpha B) (SL3-3 enhancer factor 1 a)
Human Swiss Prot No.:	Q01196
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat RUNX1 (Acetyl Lys24)
Storage and Stability:	-20°C/1 year