
Antibodies Product Datasheet

Product Name:	AR (Acetyl-Lys632) Antibody
Catalog Number:	D12224
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AR (Acetyl Lys632)
Concentration:	1 mg/ml
Gene Name:	AR DHTR NR3C4
Protein Name:	AR (Acetyl Lys632)
Other Name:	Androgen receptor (Dihydrotestosterone receptor) (Nuclear receptor subfamily 3 group C member 4)
Human Swiss Prot No.:	P10275
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse AR (Acetyl Lys632)
Storage and Stability:	-20°C/1 year