

Antibodies Product Datasheet

Product Name:	Stat3 (Acetyl-Lys49) Antibody
Catalog Number:	D12195
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Stat3 (Acetyl-Lys49)
Concentration:	1 mg/ml
Gene Name:	STAT3 APRF
Protein Name:	Stat3 (Acetyl-Lys49)
Other Name:	Signal transducer and activator of transcription 3 (Acute-phase response factor)
Human Swiss Prot No.:	P40763
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous acetyl levels of Stat3 (Acetyl-Lys49) at Human:K49, Mouse:K49, Rat:K49
Storage and Stability:	-20°C/1 year